

## PBW254-106 - September 20, 2017

Item # PBW254-106 was discontinued on September 20, 2017. For informational purposes, this is a copy of the website content at that time and is valid only for the stated product.

### CO<sub>2</sub> LASER BREWSTER POLARIZER

- ▶ 10.6  $\mu\text{m}$  Brewster Polarizer
- ▶ Zero Beam Displacement
- ▶  $\text{\O}1/2$ " and  $\text{\O}1$ " Size Options



PBW127-106

PBW254-106



PBW254-106 CO<sub>2</sub>  
Laser Polarizer  
Mounted in a KS1RS  
Kinematic Rotation Mount

#### OVERVIEW

##### Features

- Two ZnSe Plates at Brewster Angle with Polarization Coating for 10.6  $\mu\text{m}$
- Chevron Configuration Corrects for the Displacement of Transmitted Beam
- $\text{\O}1/2$ " and  $\text{\O}1$ " Versions Available

Thorlabs' CO<sub>2</sub> Laser Polarizers are designed specifically as Brewster polarizers for CO<sub>2</sub> lasers (10.6  $\mu\text{m}$ ). Each polarizer is composed of two zinc selenide windows with a polarizing coating on the outside face of each window, making this device bidirectional. The polarizer is then housed in a  $\text{\O}1/2$ " or  $\text{\O}1$ " cylindrical mount. The ZnSe windows are arranged in a chevron configuration within the cylindrical mount, resulting in no spatial offset to the transmitted beam, while the reflected s-polarization beam is ejected through slots in the cylindrical mount. The zero spatial offset allows this polarizer to be inserted and removed from a system without significantly affecting the downstream beam path. The movie to the right illustrates a beam propagating through the polarizer.

These polarizers are compatible with the SM05P05 ( $\text{\O}1/2$ " ) or the SM1P1 ( $\text{\O}1$ " ) adapters, allowing the Brewster polarizers to be mounted on any SM05- or SM1-threaded mount, respectively. Additionally, the Brewster polarizers may be mounted into a slotted lens tube (SM05L30C for  $\text{\O}1/2$ " or SM1L30C for  $\text{\O}1$ " ). The rightmost image shows the PBW254-106 in an SM1L30C slotted lens tube mounted on a KS1RS kinematic rotation mount, while the near right image shows the PBW127-106 in an SM05L30C slotted lens tube mounted on a KS05RS kinematic rotation mount. The slotted lens tube gives the option to block or unblock the beam ejected from the polarizer. However, caution should be exercised when blocking portions of the CO<sub>2</sub> beam in order to avoid potential damage to other parts.

As a result of their design, the polarizers offered here are sensitive to impact, and care should be exercised when handling these devices. Dropping a polarizer may adversely affect its performance or crack the ZnSe windows. Additionally, when handling optics, one should always wear gloves. This is especially true when working with zinc selenide, as it is a hazardous material. For your safety, please follow all proper precautions, including wearing gloves when handling these polarizers and thoroughly washing

Animation showing light propagation through Brewster Polarizer



Click to Enlarge  
PBW254-106 housed in an SM1L30C  
slotted lens tube and mounted in a  
KS1RS kinematic rotation mount



Click to Enlarge  
PBW127-106 housed in an SM05L30C  
slotted lens tube and mounted in a  
KS05RS kinematic rotation mount

your hands afterward.


Specifications	PBW127-106	PBW254-106
Design Wavelength	10.6 $\mu\text{m}$	
Clear Aperture	$\varnothing 0.14''$ (3.5 mm)	$\varnothing 0.32''$ (8.0 mm)
Polarizing Coating	$T_p \geq 95\%$ , $R_s \geq 97\%$	
Angle of Incidence <sup>a</sup>	67.4°	
Extinction Ratio <sup>b</sup>	$\geq 1500:1$	
Surface Quality	40-20 Scratch-Dig	
Surface Flatness	$\lambda @ 633 \text{ nm}$	

- The angle of incidence between laser beam and ZnSe windows. See *Graphs* tab for more information
- The extinction ratio (ER) is the ratio of maximum to minimum transmission of a linearly polarized input. When the transmission axis and input polarization are parallel, the transmission is at its maximum; rotate the polarizer by 90° for minimum transmission.

## S P E C S

Specifications	PBW127-106	PBW254-106
Design Wavelength	10.6 $\mu\text{m}$	
Polarizer Length	1.23" (31.3 mm)	2.80" (71.1 mm)
ZnSe Plate Dimension (H x W)	26 mm x 12.1 mm (1.02" x 0.48")	12 mm x 6 mm (0.47" x 0.24")
Dimensional Tolerance	$\pm 0.25 \text{ mm}$ ( $\pm 0.01''$ )	
ZnSe Plate Thickness	2.0 mm (0.08")	
Thickness Tolerance	+0.00 mm/-0.10 mm (+0.00"/-0.004")	
Polarizer Clear Aperture	$\varnothing 0.14''$ (3.5 mm)	$\varnothing 0.32''$ (8.0 mm)
Polarizing Coating	$T_p \geq 95\%$ , $R_s \geq 97\%$	
Angle of Incidence <sup>a</sup>	67.4°	
Extinction Ratio <sup>b</sup>	$\geq 1500:1$	
Surface Quality	40-20 Scratch-Dig	
Surface Flatness	$\lambda @ 633 \text{ nm}$	

- The angle of incidence between laser beam and ZnSe windows. See *Graphs* tab for more information
- The extinction ratio (ER) is the ratio of maximum to minimum transmission of a linearly polarized input. When the transmission axis and input polarization are parallel, the transmission is at its maximum; rotate the polarizer by 90° for minimum transmission.




Click to Enlarge

Click to Enlarge


## G R A P H S

### ZnSe Transmission Data




Click to Enlarge  
The transmission curve above was obtained using a 6.3 mm thick, uncoated sample of ZnSe; the incident light was normal to the surface. Please

note that this is  
the measured transmission, including surface reflections.



Click to Enlarge  
Percent  
transmission of p-  
polarized  
throughput beam.  
AOI is measured  
from the normal of  
the polarizer's face  
plate.



Click to Enlarge  
Extinction ratio of p-polarized throughput beam. AOI is  
measured from the normal of the polarizer's face plate.

Part Number	Description	Price	Availability
PBW127-106	Ø1/2" ZnSe Brewster Window Polarizer, 10.6 µm	\$1,739.00	Today
PBW254-106	Ø1" ZnSe Brewster Window Polarizer, 10.6 µm	\$1,600.00	Lead Time