


Fiber Cleavers



CAC400
Compact Cleaver



LDC401
Standard Cleaver

vytran® Fiber Cleavers

Thorlabs' lines of Vytran® Fiber Cleavers can cleave standard and specialty fibers with claddings from $\text{\O}60 \mu\text{m}$ up to $\text{\O}1.25 \text{ mm}$. Flat cleave only or flat/angled cleave units give you the flexibility to choose the best fit for your application.

These fiber cleavers feature a diamond cleave blade for industry leading longevity and precision. Each unit uses the tension-and-scribe method to achieve high-quality cleaves.

Thorlabs carries two lines of basic fiber cleavers: Standard and Compact. The Compact line of cleavers can cleave fibers with claddings from $\text{\O}60 \mu\text{m}$ to $\text{\O}600 \mu\text{m}$ while the Standard line can handle fibers with claddings from $\text{\O}80 \mu\text{m}$ up to $\text{\O}1.25 \text{ mm}$. Both lines offer flat cleave only or flat/angled cleave units.

THORLABS

Fiber Cleaving

Standard Cleavers


Our standard cleaver line can cleave fibers with claddings from $\text{\O}80\ \mu\text{m}$ to $\text{\O}1.25\ \text{mm}$. Each unit comes with a handset controller for setting cleave parameters and a built-in micrometer backstop to support specialty fibers during low-tension cleaves. These cleavers use standard Vytran fiber holder inserts, making these units easy to integrate with other Vytran workstations.

Compact Cleavers


The compact cleaver line cleaves fibers with claddings from $\text{\O}60\ \mu\text{m}$ to $\text{\O}600\ \mu\text{m}$. These cleavers are configured and controlled from a web interface that can be accessed from any web browser. An ethernet cable is included for communication.

Features

- ◆ Flat (0°) or Angled Cleaves (up to 15°)
 - ◆ Cleave Fibers of Various Sizes
 - Standard Cleaver: Claddings from $\text{\O}80\ \mu\text{m}$ to $\text{\O}1.25\ \text{mm}$
 - Compact Cleaver: Claddings from $\text{\O}60\ \mu\text{m}$ to $\text{\O}600\ \mu\text{m}$
 - ◆ Programs for Cleaving Standard or Specialty Fiber:
 - Single Mode, Multimode, and Polarization-Maintaining Fiber
 - Photonic Crystal Fiber^a
 - Capillary Tubes
 - Microstructured Fibers^a
 - Non-Circular Fibers^a
 - ◆ Mirror-Quality End Face Finishes for High-Performance Splicing
 - ◆ Robust Diamond Scribe for Longevity and High-Volume Fiber Processing
 - ◆ Industry Preferred Unit for Production Line Environments
- a. It may not be possible to create clean, angled cleaves in specialty fibers with large "air-fill" fractions due to the structural properties of the fiber.



Standard Cleaver, Flat and Angled Cleave Unit





Compact Cleaver, Flat and Angled Cleave Unit

Tension-and-Scribe Cleave Process

Thorlabs' fiber cleavers use the "tension-and-scribe" cleaving process, where tension is applied along the length of the fiber followed by an automatic scribing process using a diamond blade. Tension is maintained, causing the scribe to propagate across the fiber width and complete the cleave. Angled cleaves are accomplished by using the rotation stage to apply torsion to the fiber prior to commencing the "tension-and-scribe" process.

Certain specialty fibers, such as photonic crystal fiber, microstructured fibers, or highly stressed fibers, such as polarization-maintaining fiber, may require special parameters in order to create clean cleaves at the desired angle. For these cleaves, the initial tension applied to the fiber is lower and is slowly increased after the scribe to propagate it across the fiber and complete the cleave.




Specifications

Item #	CAC400	CAC400A	LDC401	LDC401A
Fiber Cladding	Ø60 µm to Ø600 µm		Ø80 µm to Ø1.25 mm	
Fiber Buffer/Coating	Ø67 µm to Ø1200 µm		Ø80 µm to Ø3.198 mm	
Cleave Angle	0°	0° to 15°	0°	0° to 15°
Cleave Angle Accuracy	±1.0°	±1.0°	±0.5°	±0.5° (Flat); ±1.0° (Angled)
Cleave Method	Tension and Scribe			

Standard Cleavers

Required Items for Basic Setup

- ◆ One LDC401 (Flat Cleave Only) or One LDC401A (Flat/Angled Cleave)
- ◆ Two Fiber Holder Top Inserts
- ◆ Two Fiber Holder Bottom Inserts
- ◆ Handset Controller for Large Diameter Cleavers
 - Intuitive GUI and Capacitive Touchscreen



Our digital microscope kits allow the user to image the fiber, cleave blade, and micrometer position during the cleaving process.

Additional Customization

- ◆ Top and Bottom Inserts for Fiber Holding Blocks
- ◆ Transfer Bottom Inserts for Fiber Outer Diameters ≤1047 µm
 - Corresponding Clamp for Transferring Fiber Between Processing Stations in Assembly
 - Graphite Tips for Supporting Cladding or Buffer During Transfer
- ◆ Digital Microscope Kit
- ◆ Replacement Diamond Blade



Compact Cleavers

Required Items for Basic Setup

- ◆ One CAC400 (Flat Cleave Only) or One CAC400A (Flat/Angled Cleave)
- ◆ One V-Groove Fiber Clip
- ◆ One Top V-Groove Fiber Holder Insert
- ◆ One Bottom V-Groove Fiber Holder Insert



VYT300C Handset Controller (Included with Standard Cleavers)



Additional Customization

- ◆ Eight Sizes of V-Groove Fiber Clips
- ◆ Six Top and Bottom V-Groove Insert Options
- ◆ Micrometer Backstop and Camera Assembly
 - Supports Large-Diameter or Specialty Cleaves
 - Provides 200X Magnification of Cleave Process
- ◆ Replacement Diamond Blade

Contact Us

Contact Vytran for assistance in selecting components for your specific application.

1-732-972-2880 or
techsupport@thorlabs.com



Robert Walz
Vytran General Manager

Worldwide Support



Thorlabs, Inc.
Newton, New Jersey
Phone: 1-973-300-3000
Email: sales@thorlabs.com

Thorlabs Vytran® Division
Morganville, New Jersey
Phone: 1-973-300-3000
Email: sales@thorlabs.com

Thorlabs Measurement Systems (TMS) - NJ
Blairstown, New Jersey
Phone: 1-908-362-6200
Email: tms-sales@thorlabs.com

Thorlabs Measurement Systems (TMS) - NH
Londonderry, New Hampshire
Phone: 1-973-300-3000
Email: tms-sales@thorlabs.com

Thorlabs Quantum Electronics (TQE)
Jessup, Maryland
Phone: 1-973-300-3000
Email: sales-TQE@thorlabs.com

Thorlabs Imaging Systems
Sterling, Virginia
Phone: 1-703-651-1700
Email: imagingsales@thorlabs.com

Thorlabs Spectral Works (TSW)
West Columbia, South Carolina
Phone: 1-973-300-3000
Email: sales@thorlabs.com

Thorlabs Ultrafast Optoelectronics
Ann Arbor, Michigan
Phone: 1-973-300-3000
Email: sales@thorlabs.com

Thorlabs Laser Division - CO
Lafayette, Colorado
Phone: 1-973-300-3000
Email: sales@thorlabs.com

Thorlabs Crystalline Solutions (TCS)
Santa Barbara, California
Phone: 1-973-300-3000
Email: sales@thorlabs.com

Thorlabs Canada
Phone: 1-973-300-3000
Email: sales@thorlabs.com

Thorlabs Ltda, Brazil
Phone: +55 (16) 3413 7062
Email: brasil@thorlabs.com

Thorlabs Ltd.
Phone: +44 (0)1353 654440
Email: sales.uk@thorlabs.com

Thorlabs SAS France
Phone: +33 (0) 970 444 844
Email: sales.fr@thorlabs.com

Thorlabs GmbH / Thorlabs Lübeck
Phone: +49 (0) 8131 5956-0
Email: europe@thorlabs.com

Thorlabs Elliptec® GmbH
Phone: +44 (0)1353 654440
Email: sales.de@thorlabs.com

Thorlabs Vytran® Europe
Phone: +44 (0) 1392-445777
Email: vytran.uk@thorlabs.com

Thorlabs Sweden AB
Phone: +46 31 733 30 00
Email: scandinavia@thorlabs.com

Thorlabs China Ltd.
Phone: +86 (0)21-60561122
Email: chinasales@thorlabs.com

Thorlabs Japan
Phone: +81-3-6915-7701
Email: sales@thorlabs.jp