


APPLICATION



Xenopus Coronal Plane: Ceratohyal Cartilage^{1,*}

Animal models are studied to understand biological phenomena and transfer the findings to human biology and medicine.


The non-invasive nature of OCT has made it an indispensable tool that allows researchers to image animal models *in-vivo* over the course of the animal's life into adulthood.

QUICK FACTS

- ◆ Animal models can be imaged *in-vivo* and at various stages in their lives.
- ◆ Long wavelengths such as 1300 nm penetrate deep into tissue.
- ◆ Shorter wavelengths such as 880 nm allow high-resolution imaging.
- ◆ The penetration depth may vary depending on the tissue type.
- ◆ M-modes (depth scans vs. time) highlight changes at a specific lateral position.

TYPICAL SETUP

For *in-vivo* imaging, zebrafish and xenopus are typically anesthetized and then immobilized, e.g. in agarose. OCT experiments can then be performed from different angles.^{1,2}



Protocol for xenopus immobilization and imaging.^{1,*}

After anesthetization, drosophila flies can be fixed using adhesive. OCT imaging can be performed after the flies wake up.⁵


To image chicken embryos, a window is cut into the eggshell and the vitelline membrane is peeled away. Inserting a glass window preserves the embryo and allows long-term measurements.⁶

PUBLICATIONS


- 1) E. Deniz, S. Jonas, M. Hooper, J.N. Griffin, M.A. Choma, M.K. Khokha, *Sci. Rep.*, **7**, 42506, 2017
- 2) P. Date, P. Ackermann, C. Furey, I.B. Fink, S. Jonas, M.K. Khokha, K.T. Kahle, E. Deniz, *Sci. Rep.*, **9**, 6196, 2019
- 3) Y.S. Lin, C.C. Chu, P.H. Tsui, C.C. Chang, *J. Biophotonics*, **6** (9), 668, 2013
- 4) J. Zhang, W. Ge, Z. Yuan, *Biomed. Opt. Express*, **6** (10), 3932, 2015
- 5) A. Lam, P. Karekar, K. Shah, G. Hariharan, M. Fleishman, H. Kaur, H. Singh, S. Gururaja Rao, *Sci. Rep.*, **8**, 6910, 2018
- 6) K. Courchaine, M.J. Gray, K. Beel, K. Thornburg, S. Rugonyi, *J. Cardiovasc. Dev. Dis.*, **6** (1), 11, 2019
- 7) M. Marrese, N. Antonovaite, B.K.A. Nelemans, T.H. Smit, D. Iannuzzi, *Acta Biomater.*, **97**, 524 2019 (in press)

* Images licensed under a Creative Commons Attribution 4.0 International License, see <http://creativecommons.org/licenses/by/4.0/>. Image adapted from Original.


EXAMPLE IMAGES



M-scan of drosophila heartbeat. Highlighted regions show the end diastolic diameter (EDD), end systolic diameter (ESD), diastolic interval (DI), and systolic interval (SI).^{5,*}




Top: Ventral Three Chamber View during Xenopus Cardiac Cycle^{1,*}
Right: Morphology of Chicken Embryo (3D View and Sagittal as Well as Transverse Cross Sections)^{7,*}



RECOMMENDED ITEMS

Choice of OCT System:

- ◆ **TEL221C1**: For Deep Penetration
- ◆ **GAN332C1**: For High Axial & Lateral Resolution
- ◆ **GAN632C1**: For High Speed and High Axial & Lateral Resolution



Useful Accessories:

- ◆ Mitutoyo Lenses for Very High Lateral Resolution:
 - ◆ **OCT-LKM10-SP1** & **OCT-RAM20-SP1**
4 µm @ 900 nm and 6 µm @ 1300 nm
 - ◆ **OCT-LKM20-SP1** & **OCT-RAM20-SP1**
2 µm @ 900 nm and 3 µm @ 1300 nm
- ◆ Higher Lateral Resolution Lenses Available on Request

Interested? Email OCT@thorlabs.com for more information.