

CHAPTERS

Tables/
Breadboards

Mechanics

Optomechanic
Devices

Kits

Lab Supplies

SECTIONS

Mini Series

Ø1/2" Post
AssembliesØ1" Post
AssembliesØ1.5" Post
AssembliesMounting/Angle
Brackets

Lab Platforms

Instrument Shelves

Lens Tubes

Cage Systems

Optical Rails

Mirror Mounts

Kinematic Mounts

Fixed Mounts

Translation Mounts

Rotation Mounts


Filter Mounts

V-Mounts

Iris Apertures

Adapters

Stepper and DC Servo Kits: 1/4"-80 Actuator Mirror Mounts



Combine our TST001 T-Cube Stepper Motor Driver with either our ZST6 or ZST12 Stepper Motor Actuators to motorize any mirror mount with 1/4"-80 actuators.



TPS008




TPS001




TCH002


Choose from three power supply options to upgrade your mirror mount to a motorized version.




Z806



ZST6



Z812



ZST13

DC Servo Actuators

ITEM #	\$	£	€	RMB	BACK PLATE THICKNESS	DESCRIPTION
Z806*	\$ 475.00	£ 342.00	€ 413,25	¥ 3,785.75	13 - 15 mm	6 mm Travel Motorized DC Servo Actuator, 1/4"-80 Thread
Z812**	\$ 485.00	£ 349.20	€ 421,95	¥ 3,865.45	15 - 20.5 mm	12 mm Travel Motorized DC Servo Actuator, 1/4"-80 Thread

*Recommended for use with our KM100 2-Adjuster Kinematic Mirror Mount

**Recommended for use with our KS1 3-Adjuster Kinematic Mirror Mount

Stepper Motor Actuators

ITEM #	\$	£	€	RMB	BACK PLATE THICKNESS	DESCRIPTION
ZST6*	\$ 576.30	£ 414.94	€ 501,38	¥ 4,593.11	13 - 15 mm	6 mm Travel Motorized Stepper Motor Actuator, 1/4"-80 Thread
ZST13**	\$ 601.80	£ 433.30	€ 523,57	¥ 4,796.35	15 - 20.5 mm	13 mm Travel Motorized Stepper Motor Actuator, 1/4"-80 Thread

*Recommended for use with our KM100 2-Adjuster Kinematic Mirror Mount

**Recommended for use with our KS1 3-Adjuster Kinematic Mirror Mount

DC Servo and Stepper Motor Controllers (Power Supply Sold Separately)

ITEM #	\$	£	€	RMB	DESCRIPTION
TDC001	\$ 595.00	£ 428.40	€ 517,65	¥ 4,742.15	T-Cube DC Motor Controller
TST001	\$ 595.00	£ 428.40	€ 517,65	¥ 4,742.15	T-Cube Stepper Motor Controller

Power Supply Options

ITEM #	\$	£	€	RMB	DESCRIPTION
TPS001	\$ 25.00	£ 18.00	€ 21,75	¥ 199.25	15 V Power Supply for a Single T-Cube Unit
TPS008	\$ 175.00	£ 126.00	€ 152,25	¥ 1,394.75	15 V Power Supply for up to Eight T-Cube Units
TCH002	\$ 726.90	£ 523.37	€ 632,40	¥ 5,793.39	T-Cube Controller Hub and Power Supply Unit



Combine our TDC001 T-Cube DC Motor Driver with either our Z806 or Z812 DC Motor Actuators to motorize any mirror mount with 1/4"-80 actuators.