

## 19" Equipment Racks


### Standard 19" EIA Rack

- 7' x 19" Floor-Standing Equipment Rack
- Heavy-Duty, 3" Wide, Vertical Channels
- Weight Capacity is 750 lbs (Load Must be Evenly Distributed)
- 12-24 Tapped Mounting Holes in a 5/8", 5/8", 1/2" Repeating Hole Pattern
- Side Holes Allow Multiple Racks to be Bolted Together
- Universal Design is Imperial and Metric Compatible

RK3884



These heavy-duty, 19" equipment racks are an ideal lab organizer. Manufactured from aluminum U-channels and angle extrusions, the weight capacity is 750 lbs, which should be distributed evenly along its height.



ITEM #	\$	£	€	RMB	DESCRIPTION
RK3884	\$ 248.00	£ 178.56	€ 215.76	¥ 1,976.56	7' x 19" Wide Equipment Rack

## 19" Rack Mounting Shelves



RK4101

Shown with Essentials Kit  
(See Pages 374 - 380 for Details)




RK4102

The RK4102 has 77 Ø1/2" holes on a 1.5" grid. This allows components on Ø1/2" post (see pages 93 - 94) to be conveniently stored.

RK3884

19" Instrumentation Rack



RK4101

The RK4101 is a great way to organize our Essentials Kits within a standard 19" rack. Once mounted, these shelves provide a convenient way to store optical elements.



ITEM #	\$	£	€	RMB	DESCRIPTION
RK4101	\$ 59.00	£ 42.48	€ 51.33	¥ 470.23	Plain Fixed Rack Shelf, 17.25" x 11.20" Surface
RK4102	\$ 84.00	£ 60.48	€ 73.08	¥ 669.48	Fixed Rack Shelf, Ø1/2" Holes, 1.5" Centers