

1/4"-20 Coarse Hex Adjusters


Features

- 1" and 2" Screw Lengths
- 0.57" Long Bushings
- Screws Have a Hex Head (5/64" Hex Key)
- Removable Knurled Knob for Easy Adjustment
- Vacuum-Compatible Lock Nut

These 1/4"-20 thread adjustment screws are for use in assemblies where fine adjustment is not necessary. They are supplied in 1" and 2" lengths. Bushings, lock nut, and a knurled adjustment knob are also available.



Mechanical
Drawings Available on the
WEB



Please refer to our website for complete models and drawings.

ITEM#	\$	£	€	RMB	DESCRIPTION	LENGTH
CS25SS100	\$ 6.20	£ 4.30	€ 5.60	¥ 52.40	1/4"-20 Coarse Adjustment Screw	1"
CS25SS200	\$ 9.20	£ 6.40	€ 8.20	¥ 77.70	1/4"-20 Coarse Adjustment Screw	2"
CS25SSN1	\$ 5.70	£ 3.95	€ 5.10	¥ 48.20	1/4"-20 Bushing	0.57"
CS25SSK1	\$ 7.20	£ 5.00	€ 6.40	¥ 60.80	1/4"-20 Removable Adjustable Knob	0.47"
LN2520	\$ 8.20	£ 5.70	€ 7.30	¥ 69.30	1/4"-20 Lock Nut	-

Polaris™ Kinematic Mirror Mount

Pointing the Way to Precision Alignment

Design Elements

- Thermal and Dimensional Stability: Heat-Treated Stainless Steel
- Repeatability and Durability: Ball Contacts with Sapphire Seats
- Usability and Feel: Matched Actuator/Bushing Pairs

Temperature Cycle Testing:
Less than 2 μrad Final Deviation
Additional Mechanical and Temperature Test Data at www.Thorlabs.com

NEW!



See Page 214 for More Details