

## Mini-Series Table Clamps



Our MSCL series of mini breadboard clamps have been fashioned in the same spirit as our popular CL series clamps (see page 82). These versatile clamps provide excellent flexibility and secure clamping force for positioning components at random across the mini breadboard featured on page 17.



ITEM#	METRIC ITEM#	\$	£	€	RMB	DESCRIPTION
MSCL1	MSCL1*	\$ 6.10	£ 3.80	€ 5,70	¥ 58.30	0.8" (20mm) Mini Table Clamp
MSCL2	MSCL2/M	\$ 7.10	£ 4.50	€ 6,60	¥ 67.80	1.6" (40mm) Mini Table Clamp
MSCL3	MSCL3*	\$ 7.10	£ 4.50	€ 6,60	¥ 67.80	1.2" (29mm) Mini Table Clamp

\*Imperial and Metric Compatible

## Mini-Series Post Clamps

The Mini-Series MSRA90 provides 90° alignment between two 6mm Mini-Series Rods. The through-holes are double-bored for two lines of post-to-post contact. Thumbscrews located on the top and bottom of this adapter secure the rods in place.

- Compatible With 6mm Mini Series Rods
- Provides Control and Flexibility in Tight Setups
- Essential Building Blocks for Complex, Rigid, 3D Structures

The ER90B right angle post adapter is used to mount two of our new Mini-Series posts at right angles to each other. The adapter is double-bored to provide two stable lines of contact between the mating post and adapter body. This simple adapter is an essential building block for constructing optical systems from our new family of miniature components.

The MSWC Mini-Series swiveling post clamp provides 360° of continuous adjustment. This feature allows any arbitrary angle to be achieved, making even the most complicated construction simple. Once the desired angle has been achieved, a locking screw secures the swiveling mechanism. Each hole is double-bored for two lines of stable contact.



ITEM#	METRIC ITEM#	\$	£	€	RMB	DESCRIPTION
MSRA90	MSRA90*	\$ 15.00	£ 9.50	€ 14,00	¥ 143.30	Right Angle Post Clamp - Mini Series
ER90B	ER90B/M	\$ 12.40	£ 7.80	€ 11,50	¥ 118.40	90° Mini Series Post Adapter
MSWC	MSWC*	\$ 33.00	£ 20.80	€ 30,70	¥ 315.20	90° Mini Series Swiveling Post Clamp

\*Imperial and Metric Compatible

## Ø6mm Post Thread Adapter

Our MSA8 and MSA25 adapters can be used to adapt Ø6mm posts to other thread standards. Machined from cold rolled steel, the adapters can be used with our Mini-Series Posts (see page 76).



ITEM#	METRIC ITEM#	\$	£	€	RMB	DESCRIPTION
MSA8	MSA4/M	\$ 7.99	£ 5.00	€ 7,40	¥ 76.30	6mm Post to #8-32 (M4) Adapter
MSA25	MSA6/M	\$ 7.99	£ 5.00	€ 7,40	¥ 76.30	6mm Post to 1/4-20 (M6) Adapter

\*Imperial and Metric Compatible

## Ø6mm Post Mini-Series Collar

The MSR2 6mm slip-on post collar is a mini version of our popular R2 Ø1/2" post collars. It features a double-bore profile to provide a stable two-line contact, and a non-marring brass setscrew to ensure secure locking. Height adjustment is restricted while providing rotational freedom.



MSR2 SLIP-ON COLLAR  
Mini-Series Post Sold Separately

ITEM#	\$	£	€	RMB	DESCRIPTION
MSR2	\$ 6.00	£ 3.80	€ 5,60	¥ 57.30	Mini-Series Post Collar