

Lab Supplies

Lens Care/Cleaning

Epoxies

Cables, Electrical


Safety & Blackout

Screws

19" Rack Hardware

Tools

Optic Storage Boxes



Sold in Packages of 10

ITEM#	\$*	£*	€*	RMB*	DESCRIPTION
BX01	\$ 16.00	£ 10.10	€ 14.90	¥ 152.80	Optic Box for Ø1" Optics (Pack of 10)
BX02	\$ 21.00	£ 13.20	€ 19.50	¥ 200.60	Optic Box for Ø2" Optics (Pack of 10)
BX03	\$ 14.00	£ 8.80	€ 13.00	¥ 133.70	2.5" x 2.5" x 1" Boxes w/ Foam Insert

*Price is per package of 10.

Inspection Loupe and Magnifier



JEL10

- 10X Magnification
- Hastings Design
- Lightweight

Jeweler's Eye Loupe

This Bausch & Lomb loupe is ideal for inspecting optics and small parts. It provides a clear image and allows for hands-free examination.

Head-Worn Magnifier

This head-worn inspection magnifier is ESD compliant and cleanroom compatible. It is ideal for inspecting components when high magnification is not required, and may be worn with or without eyeglasses. The magnifying lenses can be swung out of the way or locked into position.



MAG200K

ITEM#	\$	£	€	RMB	DESCRIPTION
JEL10	\$ 15.86	£ 10.00	€ 14.70	¥ 151.50	Eye Loupe, 10X
MAG200K	\$ 29.70	£ 18.70	€ 27.60	¥ 283.60	ESD Inspection Magnifier, 2.0X
MAG275L	\$ 13.95	£ 8.80	€ 13.00	¥ 133.20	Replacement Lens for MAG200K, 2.75X


- 2.0X (2.75X With Replacement Lenses) Magnification
- ESD Compliant
- Worn With or Without Eyeglasses

Inspection Mirrors



TM1

- Corrosion-Resistant Stainless Steel Telescoping Handle
- High-Quality, Distortionless, Glass Mirror
- Articulated Joint Adds Versatility in Mirror Positioning
- 8.50" Solid Length, 35" Extended Length



Partial Extension of Telescoping Arm

ITEM#	\$	£	€	RMB	DESCRIPTION
TM1	\$ 20.00	£ 12.60	€ 18.60	¥ 191.00	Inspection Mirrors

Multilayered Adhesive Mats

Features:

- 18" x 36"
- 30 Sheets
- Acrylic Adhesive Polyethylene Sheets



Multilayered Adhesive Mats effectively capture dirt and dust from foot traffic and equipment wheels before entering a controlled environment. Each layer is numbered, which helps ensure one sheet removal and identify the number of remaining applications. Low-profile construction reduces the risk of tripping and allows free movement of foot or wheel traffic.



ESD20

ITEM#	\$	£	€	RMB	DESCRIPTION
ESD20	\$ 52.00	£ 32.80	€ 48.40	¥ 496.60	Adhesive Floor Mat